

# ODEL

ENHANCING **LIFE**  
THROUGH **TECHNOLOGY**

## Environmental Control

ODEL are able to provide a wide range of environmental control devices to meet all needs. These can range from simple 4 button devices to complex devices to operate many home appliances and integrate with wheelchair controls.

### HouseMate Life

HouseMate Lite is a stand-alone switch accessible infra-red remote control with auditory feedback using recorded speech. Practical and easy to use, HouseMate can be operated by external switches or directly by pressing the keypad.



### HouseMate Pro

The HouseMate Pro, allows multifunction access to an Android mobile phone by the simple clicking of an external Bluetooth switch, joystick or by using the touchscreen of the Android mobile. Any single switch or 5 switch device can be used as well as a powered wheelchair joystick. Virtually all of the phone's functionality can be accessed through the software, making and answering calls, sending and replying to texts, browsing the Internet, send and receive e-mails, configure telephone settings, access the weather and markets as well as access to most available apps.

info@odel.co.uk  
Banbridge: 028 4062 8349

Buy online at...  
[www.odelmobility.co.uk](http://www.odelmobility.co.uk)



## The Relax

### Accessible IR Learning Remote

The NEW Relax is an accessible infrared (IR) learning remote that comes with an updated modern design and provides simplified control of almost any IR device. With a light touch or with the use of one switch, an individual can control up to eight functions on their favourite IR device.

#### Key Features:

- Switch Accessible
- Simplified Design
- Programmable
- Mountable



## PICOMED Devices

The PICO 4/8/16 are IR transmitters with 2 levels and with scanning. They have built in IR-signals and are also able to learn other IR-signals. They come with a 3.5mm mono jack socket, to enable connection of a switch.

With ergonomically designed, soft keys with clearly visible icons and built in finger guide.



## iPad Mounting Solutions

ODEL have a range of mounting solutions for iPad devices. These range from desktop stands, wheelchair mounts and many more. For more information on our full range, visit us at [www.odelmobility.co.uk](http://www.odelmobility.co.uk)



## Switch Solutions

ODEL can provide a wide range of switch solutions to meet all needs. Here is just a small selection.



Buy online at...  
[www.odelmobility.co.uk](http://www.odelmobility.co.uk)